

B2M Control Normal

Volume: 4*1 ml

Material: Liquid

Preface:

B2M Calibrator is a liquid Calibrator which is prepared from human serum.

Preparation:

Open the vials lid Carefully and be aware not to lose the contents of the vial.

The control is ready to use.

Avoid shaking the vial violently and creating foam.

Storage and Stability:

All the components of vial are stable until the expiration date on the label at 2-8 °C.

Warning and precautions:

This product is only for laboratory use.

Add: NO.68, Sorena St., Damavand 2 Industrial Park, Damavand, Tehran, IRAN

Tell: 021-77778007-8

website: www.pts-ico.ir

Fax: 021-26142195

email: pts.ico@gmail.com

V.F1